



- compact active electrode for high-density EEG and multi-purpose application
- fully compatible with g.GAMMAsys and g.Hlamp's active headbox
- touchproof sintered Ag/AgCl-ring electrode in a durable molded housing
- super-fast and reliable setup: fill-and-ready
- usable with g.GAMMAcap² or with adhesive washers
- simultaneous impedance check of all channels with g.Hlamp





Electrode holder rings with channel labels attached to g.GAMMAcap²

Adhesive washers can be used to mount the electrode holder ring directly on the skin (EOG, ECG, EMG, ...)



g.GAMMAgel or g.GAMMAcream

can be used to fill the electrodes with a syringe equipped with a 2-mm blunt cannula (no scratching or skin abrasion required)



No waste of time and no mistakes: electrodes can remain in the cap and plugged in the connector box for cleaning

> Easy cleaning with a special brush



The g.SCARABEO electrode snaps into the electrode holder ring, labels can be used to identify channels and positions

> A durable molded body protects the touchproof Ag/AgCl ring electrode and seals the sensitive electronics





A small reservoir for gel inside the electrode allows for many hours of recording



contact: g.tec medical engineering GmbH, Sierningstr. 14, 4521 Schiedlberg, AUSTRIA phone: +43 7251 22240, fax: +43 7251 22240 39, e-mail: office@gtec.at, web: www.gtec.at